DTVNotification

To: Director, Clinical/Biomedical Engineering Engineer, Risk Manager, Hospital

Administrator, Nursing Home Director

From: WNYW New York, NY Fox Television Station, LLC

Date: July 3rd 2018

Notice of Facility Change for Digital Television Station

You are receiving this letter because WNYW in New York, New York has received authority from the Federal Communications Commission and will be changing the location, antenna pattern, and power of its transmitter which operates on UHF Channel 44 (650-656 MHz). This change in transmitter technical parameters may affect equipment in your facility. Please be sure that anyone in your facility who is responsible for the supervision or maintenance of the biomedical equipment receives a copy of this notification and the attachments.

Under authorization from the Federal Communications Commission (FCC), the station listed above will be broadcasting on the frequencies with the technical parameters listed in the following pages.

Per the rules of the FCC (47 CFR 15.242) full power stations are required to make a good faith effort to identify and notify health care facilities (e.g., hospitals, nursing homes, see 47 CFR 15.242(a)(l)) within the service area potentially affected by the DTV operations.

DTV Notification has been retained by the station listed above to make the required good faith notifications to applicable facilities in their service areas. This letter is to notify health care facilities in the coverage area in the event that any electronic medical devices are affected by the radiofrequency (RF) field emitted by the transmitted DTV signals.

Your facility may have electronic biomedical devices that could be impacted by the DTV operations listed on the following pages. The FCC allows certain electronic biomedical devices to use unoccupied DTV frequencies on a secondary use basis. As a result of the secondary use status, your facility is responsible for identifying and resolving any potential radio frequency interference for the operations identified on the following pages.

To address potential interference problems for medical telemetry systems, the FDA recommends that you work with the manufacturer of your telemetry systems to:

- Determine the channel/frequencies your telemetry systems use; and
- Consult the FCC table of allotments for DTV channels.

If your telemetry systems are transmitting on channels that are scheduled to be used by a local station, (as identified in the attached sheets) you should work with the manufacturer of your telemetry systems to:

• Change your telemetry channels to unused channels prior to the expected broadcast date to avoid interference

It is recommended that you periodically check the FCC allotments to assure that those channels remain available.

https://www.fda.gov/MedicalDevices/DigitalHealth/WirelessMedicalDevices/default.htm

Sincerely,

Support Team DTVNotification.com

305 B Centennial Street, La Plata, Maryland 20646 (888)-244-0436 Ext. 1 www.dtvnotification.com

IMPORTANT MEDICAL TELEMETRY INFORMATION NOTICE

It is possible that certain medical devices in your facility may be affected by the operations identified herein. Please note, WNYW in New York, New York will commence transmitting from its new location with new technical parameters on or about August 6, 2018 on UHF Channel 44 (650-656 MHz) from a new antenna location and antenna height. With this letter we are informing you of the intent of the station to begin broadcasting on Channel 44 as authorized by the FCC with the new technical facility parameters.

As the guidance from the FDA notes:

Many medical telemetry devices are operated under 47 CFR Part 15 of the FCC regulations. These devices are secondary users of the radio-frequency (RF) spectrum. Licensed users such as television stations are the primary users. As a secondary user, your telemetry devices may be subject to interference from the primary user at any time, but must not cause interference with the primary user. If your telemetry devices are operating on frequencies licensed to a primary user, you need to be aware of the potential for interference and take any steps necessary to avoid device malfunctions due to interference

FDA Recommendations

To address potential interference problems for medical telemetry systems, FDA recommends that you work with the manufacturer of your telemetry systems to:

- Determine the channel /frequencies your telemetry systems use; and
- Consult the FCC table of allotments for DTV channels in your area.

If your telemetry systems are transmitting on channels that are or will be used by a local station, you should work with the manufacturer of your telemetry systems to:

• Change your telemetry channels to unused channels prior to the expected broadcast date to avoid interference

It is suggested that you should periodically check the FCC allotments to assure that those channel s remain available.

https://www.fda.gov/MedicalDevices/DigitalHealth/WirelessMedicalDevices/default.htm

STATION INFORMATION:

WNYW in New York, New York FCC File Number: 0000041037

FCC Facility ID: 22206

DTV channel: 44

Frequency band: 650-656 MHz

Estimated broadcast date: August 6, 2018 Antenna Effective radiated power: 151 KW

Antenna Height: 530.4 meters above mean sea level (AMSL) Antenna Location: 40° 42' 46.8"N 074° 00' 47.3 W (NAD83)

Technical Contact: Ed Harris, WNYW Engineering, Telephone: (212) 452-3638

Email: Edward.Harris@foxtv.com

DTV Notification does not represent or warrant any of the information contained in this letter. All of the information in this letter has been provided and approved by the appropriate station contact and not DTV Notification. DTV Notification further recommends that you contact the appropriate station contact listed above to confirm the information in this letter to avoid any confusion. DTV Notification shall not be held liable for any loss as a result of the information contained in this letter.

Eric Wheeler 07/03/2018



TITLE

FILE NAME

DOCUMENT ID

STATUS

WNYW Medical Letter Confirmation

WNYW Medical Lett...EDITS_070318.docx

68344543bad394f62fdf8e37a263500a9075aa8a

Completed

Document History

SENT

07/03/2018

20:09:36 UTC

Sent for signature to Rick Wheeler (Rick.Wheeler@foxtv.com)

from info@dtvnotification.com

IP: 47.215.185.6

 \odot

07/03/2018

Viewed by Rick Wheeler (rick.wheeler@foxtv.com)

VIEWED 20:34:28 UTC

IP: 8.37.180.2

r

07/03/2018

Signed by Rick Wheeler (rick.wheeler@foxtv.com)

SIGNED

20:36:37 UTC

IP: 8.37.180.2

Ø

07/03/2018

COMPLETED 20:36:37 UTC

The document has been completed.